
Technique for increasing the strength of fixation to the suspensory button in anterior cruciate ligament reconstruction using the quadriceps tendon

Department of Orthopaedic Surgery, Graduate School of Medical Sciences, Kanazawa University

Tomoyuki Kanayama, Junsuke Nakase, Yoshihiro Ishida,

Yusuke Yanatori, Yu Arima, Hiroyuki Tsuchiya

Disclosure

We have nothing to disclose

Introduction

- ✓ Quadriceps tendon (QT) has recently gained great interest
as an autograft for ACL reconstruction.
- ✓ QT exhibits better microstructural and mechanical properties
than the Hamstring tendon and Bone patellar tendon bone¹⁾.
- ✓ Primary ACL reconstruction using QT autografts
appears to have successful outcomes with a low rate of graft failure²⁾.

Fixation strategy of graft

✓ Insufficient strength for graft fixation leads to graft failure³⁾.

✓ QT cannot be folded and there is no consensus on how to fixation them.




Fixation with hamstring tendon





Fixation with QT

➤ Surgeons appear to use various fixation techniques in clinical settings.



Fixation strategy of graft

- ✓ Evaluated the elongation and load to failure of different stitching methods and suture diameters⁶.
- ✓ Direct stitch to the suspensory button device may be better from a biomechanical perspective⁷.





Purpose

To compare the biomechanical strength of various QT fixation using adjustable loops in soft-tissue QT grafts for ACL reconstruction

Material & methods


- ✓ Thirty fresh-frozen bovine Achilles tendons (n=10 in each group)
- ✓ Grafts: 10 mm width, 50 mm length, and 4 mm thick
- ✓ ACL TightRope II implant, #2 FiberWire, and FiberTag (Arthrex) were used




All knots were tied five times and fixed.

Fixation configurations


Group A



Group C



Tensile testing



AG-X plus (Shimadzu Autograph)

- ✓ **Preloading:**

 - 50 N × 5 cycles

 - statically held at 50 N for 1 min.

- **The amount of elongation (mm)**


- ✓ **Load-to-failure testing:**

 - conducted until rupture at 5mm / min.


- **Maximum load to failure (N)**

Results

Graft elongation



Maximum load to failure





To prevent graft failure


It is particularly important that the tendon

- ✓ Not damaged
- ✓ Does not stretch
- ✓ Has sufficient strength

Present study




Graft elongation



10.3 ± 2.4 mm



Group A

**Fixation on
a surface**



10.0 ± 1.0 mm

Group C




16.6 ± 2.2 mm

Group B

**Fixation in
one direction only**

Conclusion



Fixation using Speed Whip Rip Stop technique may be better from a biomechanical perspective at time zero.

Reference

- 1) Castile RM, Jenkins MJ, Lake SP, Brophy RH (2020) Microstructural and mechanical properties of grafts commonly used for cruciate ligament reconstruction. *J Bone Joint Surg* 102(22):1948–1955
- 2) Kanakamedala AC, de Sa D, Obioha OA, Arakgi ME, Schmidt PB, Lesniak BP et al (2019) No difference between full thickness and partial thickness quadriceps tendon autografts in anterior cruciate ligament reconstruction: a systematic review. *Knee Surg Sports Traumatol Arthrosc* 27(1):105–116
- 3) Hong CK, Lin CL, Chang CH, Jou LM, Su WR (2014) Effect of the number of suture throws on the biomechanical characteristics of the suture-tendon construct. *Arthroscopy* 30:1609-1615
- 4) Hughes JD, Vaswani R, Gibbs CM, Tisherman RT, Musahl V (2020) Anterior cruciate ligament reconstruction with a partial-thickness quadriceps tendon graft secured with a continuous-loop fixation device. *Arthrosc Tech* 9(5):e603–e609
- 5) Todor A, Caterev S, Nistor DV, Khallouki Y (2016) Free bone plug quadriceps tendon harvest and suspensory button attachment for anterior cruciate ligament reconstruction. *Arthrosc Tech* 5(3):e541–e544
- 6) Michel PA, Domnick C, Raschke MJ, Kittl C, Glasbrenner J, Deitermann L et al (2019) Soft tissue fixation strategies of human quadriceps tendon grafts: a biomechanical study. *Arthroscopy* 35(11):3069–3076
- 7) Kamada K, Nagai K, Nagamune K, Hoshino Y, Nakanishi Y, Araki D et al (2022) Direct suturing quadriceps tendon to a continuous loop with a suspensory button provides biomechanically superior fixation in ACL reconstruction. *Knee Surg Sports Traumatol Arthrosc* 30:2307-2313