

# Operator variation in applying a knee loading device for evaluation of tibial component loosening in total knee arthroplasty

George S. Buijs<sup>1</sup>, Maaïke A. ter Wee<sup>2</sup>, Chiel Klein<sup>1</sup>, Lidwine B. Mokka<sup>3</sup>, Johannes G. G. Dobbe<sup>2</sup>, Mario Maas<sup>4</sup>, Matthias U. Schafroth<sup>1</sup>, Geert J. Streekstra<sup>2</sup>, Leendert Blankevoort<sup>1</sup> and Arthur J. Kievit<sup>1</sup>

<sup>1</sup>Department of Orthopedic Surgery and Sports Medicine, Amsterdam UMC, Amsterdam Movement Sciences, Universiteit van Amsterdam, The Netherlands

<sup>2</sup>Department of Biomedical Engineering & Physics, Amsterdam UMC, Amsterdam Movement Sciences, Universiteit van Amsterdam, The Netherlands

<sup>3</sup>Department of Epidemiology and Data Science, Amsterdam UMC, Amsterdam Public Health, Vrije Universiteit Amsterdam, The Netherlands

<sup>4</sup>Department of Radiology & Nuclear Medicine, Amsterdam UMC, Amsterdam Movement Sciences, Universiteit van Amsterdam, The Netherlands



**ISAKOS**  
CONGRESS  
2025



**MUNICH**  
**GERMANY**  
June 8-11



**ORTHOPEDIC SURGERY  
AND SPORTS MEDICINE**



**Amsterdam UMC**

# Faculty Disclosure Information

- Leendert Blankevoort, Arthur J. Kievit and Matthias U. Schafroth are listed as inventors on a patent for a loading device and method to be used to quantify and visualise implant displacement
- The remaining authors declare no conflict of interest
- This study was funded by an internal TKI-PPP grant from the Amsterdam University Medical Center. Allowance made available by Health-Holland, Top Sector Life Sciences & Health, to stimulate public private partnerships.



- This study was funded in part by Zimmer Biomet as part of a larger cohort study. This sub-analysis was performed independent of Zimmer Biomet.





# The problem

- Aseptic loosening of the total knee arthroplasty (TKA) is an indication for revision surgery<sup>1</sup>
- Imaging modalities that are used in the diagnostic process do not measure looseness <sup>2,3</sup>
- A CT-based method was developed to measure tibial component displacement under a valgus and varus load as applied by a loading device<sup>4</sup>
- The valgus-varus loading device is applied by an operator
- Operator variation of applying the valgus-varus loading device may affect the measurement of the displacement of the tibial component

## The question

- What is the inter-operator variation, in terms of reliability and measurement error metrics, of applying the loading device for detecting tibial TKA component displacement?



# Methods – measurement of load-induced displacement

AtMoves\*  
Knee  
Loading  
Device



valgus




varus

Applied  
load  
20 Nm




Siemens Somatom Force




CT  
valgus

CT  
varus




Implant Loosening Report



AtMoves\*  
CT analysis

# Methods – Implant displacement report



## Key variables

Screw axis rotation

Maximum total point motion (MTPM)

Mean target registration error (mTRE)





# Methods

## Patient selection

- Subgroup of clinical follow-up of a Zimmer Biomet TKA
  - ✓ Vanguard PS Open Box Porous Femoral component
  - ✓ Implantation between 2009 and 2013
  - ✓ Surgery conducted at Amsterdam UMC
- Required sample size: n=16

# Methods

## Procedures

- Two immediately consecutive CT scans with loading device
  - ✓ Two independent and blinded operators
  - ✓ Identical scan protocol
- Evaluation of the CT scans
  - ✓ One operator
  - ✓ Identical evaluation protocol
- Statistical analyses<sup>5</sup>

### Components of measurement instrument

Equipment

Load application

Collection of raw CT data

Data processing & storage

Evaluation of the scans

### Study design

CT scanner

Loading device  
**operator 1**

Loading device  
**operator 2**

Technician

Technician

Automatic

Automatic

**Operator 1\***

**Operator 1\***

\* Previous results showed no influence of the operator on the assignment of the score

# Results – Patient characteristics – n=16

<b>Age in years</b> mean (SD)	59.1	(9.5)
<b>Sex</b> n (%)		
<b>Male</b>	5	(31.2)
<b>Female</b>	11	(68.8)
<b>Time interval between arthroplasty and Scan in years</b> median (IQR)	11.9	(11.4-12.5)
<b>Charnley score</b> n (%)		
<b>0:</b>	8	(50)
<b>1:</b>	4	(25)
<b>2:</b>	4	(25)
<b>ASA score</b> n (%)		
<b>1:</b>	10	(62.5)
<b>2:</b>	6	(37.5)

<b>KOOS sub scores</b>	(median, [IQR])
<b>KOOS-Pain</b>	94.4 (83.3-100)
<b>KOOS-Other Symptoms</b>	85.7 (79.4-93.8)
<b>KOOS-ADL</b>	91.9 (86.4-96.3)
<b>KOOS-Sport/Rec</b>	45.0 (20.0-70.0)
<b>KOOS-QoL</b>	75.0 (60.9-90.6)

n	Number
ASA	American Society of Anesthesiologists
KOOS	Knee Osteoarthritis Outcome Score
ADL	Activities of Daily Living
QoL	Quality of Life



# Results – error statistics

## Intraclass Correlation Coefficient

	ICC	CI
Rotation [°]	<b>0.64</b>	0.24 – 0.86
MTPM [mm]	<b>0.81</b>	0.55 – 0.93
mTRE [mm]	<b>0.84</b>	0.61 – 0.94

Rotation Rotation about screw axis [°]  
MTPM Maximum Total Point Motion [mm]  
mTRE mean Total Registration Error [mm]  
ICC Intraclass Correlation Coefficient  
CI 95% Confidence Interval


## Error Statistics

	SEM
Rotation [°]	<b>0.14</b>
MTPM [mm]	<b>0.10</b>
mTRE [mm]	<b>0.06</b>


SEM Standard Error of Measurement

# Results – Bland-Altman plots


Screw axis rotation [°]



Maximum total point motion [mm]



Mean target registration error [mm]



○ Abberant measurement explained:  
Retrospective evaluation of CT-scan revealed incorrect placement of loading device on the knee joint during the second examination





# Discussion

- Operator variation does not affect implant displacement measurement.
  - The inter-operator difference in applying the loading device resulted in a standard error of measurement (SEM) of 0.06–0.14 mm or °, which aligns with the previously established measurement error in an experimental cadaver setting (0.07–0.13 mm or °)<sup>6</sup>.
- Apart from the identified aberrant measurement, the findings do not suggest a need for further standardization of operator training.
- A single total knee arthroplasty design was included; only tibial component motion was evaluated.
- The relatively small and uniform population may limit the generalizability of the results to a broader and more diverse patient population as seen in clinical practice.

# Conclusion

The **inter-operator reliability** of applying the loading device for detecting load-induced implant displacement is **moderate** for the rotation measure and **good** for the measures maximum total point motion (MTPM) and mean target registration error (mTRE).



ISAKOS  
CONGRESS  
2025



MUNICH  
GERMANY  
June 8-11



ORTHOPEDIC SURGERY  
AND SPORTS MEDICINE



Amsterdam UMC



# References

1. Gademán, M.G.J., Van Steenberghe, L.N., Cannegieter, S.C., Nelissen, R., Marang-Van De Mheen, P.J., 2021. Population-based 10-year cumulative revision risks after hip and knee arthroplasty for osteoarthritis to inform patients in clinical practice: a competing risk analysis from the Dutch Arthroplasty Register. *Acta Orthop* 92, 280-284.
2. Barnsley, L., Barnsley, L., 2019. Detection of aseptic loosening in total knee replacements: a systematic review and meta-analysis. *Skeletal Radiology* 48, 1565-1572.
3. Buijs, G.S., Kooijenga, A.C., Rikken, R.G.H., Schafroth, M.U., Kievit, A.J. & Blankevoort, L. 2024. MRI and SPECT/CT demonstrate, with low certainty of evidence, the highest diagnostic accuracy for aseptic knee arthroplasty loosening; a systematic comparative diagnostic test review and metaanalysis. *Knee Surg Sports Traumatol and Arthrosc* 32, 2061-2074.
4. Buijs, G.S., Kievit, A.J., Ter Wee, M.A., Magg, C., Dobbe, J.G.G., Streekstra, G.J., Schafroth, M.U., Blankevoort, L. 2025. Non-invasive quantitative assessment of induced component displacement can safely and accurately diagnose tibial component loosening in patients: A prospective diagnostic study. *Knee Surg Sports Traumatol Arthrosc* 33(1):274-285.
5. Mokkink L.B., Boers M., van der Vleuten C.P.M., Bouter L.M., Alonso J., Patrick D.L, et al. 2020. COSMIN Risk of Bias tool to assess the quality of studies on reliability or measurement error of outcome measurement instruments: a Delphi study. *BMC Med Res Methodology*. 2020;20(1): 293.
6. Kievit, A.J., Buijs, G.S., Dobbe, J.G.G., Ter Wee, M.A., Kerkhoffs, G.M.M.J., Streekstra, G.J., Schafroth, M.U., Blankevoort, L. 2023. Promising results of a non-invasive measurement of knee implant loosening using a loading device, CT-scans and 3D image analysis. *Clin Biomech (Bristol)* 104:105930.